



UBILITY™

FOR PATIENTS & CAREGIVERS

MEET UbiFlex™ FROM UBILITY™


A more personal approach to arm support.

UbiFlex™ provides powered elbow assistance for patients with arm paralysis.

Brace weight: 1 kg (excluding custom-fabricated cuffs).



UbiFlex™ powered arm brace.



Muscle activity is detected as a signal to request brace activation.

Built around you.

Wheelchair use or complete paralysis of the affected arm does not rule out evaluation. A functional muscle elsewhere may provide the sEMG signal to request movement. Your clinical team assesses the signal, fit, and movement goals.

01

Talk about goals

Start with your needs and a clinical assessment.

02

Capture your fit

Use scans and measurements of your arm.

03

Create custom parts

Design and make the interfaces around you.

04

Fit and explain

Review alignment, settings, and instruction.

Questions worth asking

- Am I a candidate for the UbiFlex™?
- What would fitting and learning to use the brace involve?
- Would I need help from a caregiver?

Questions? Connect with UBILITY™.

Product support: support@ubility.inc · 330-356-5052

General inquiries: info@ubility.inc

4305 Mt Pleasant St NW, Suite 203B · North Canton, Ohio 44720

For questions about clinical suitability, speak with your orthotist or clinical team.